

# Brightree – Diagnostic Testing Integration

Oximetry & Home Sleep Testing

## Brightree Settings Setup

---

1. System Setup > User Group > Permissions > System Setup and change the Diagnostic Testing (Setup) to Full Control

\*Complete so all users can access this feature if possible.

2. System Setup > Setup and Diagnostic Testing to display the Diagnostic Test page.

3. System Setup > User Groups > Security Tab > Functional Security > Check the box for "Order Diagnostic Test"

4. If using Document Management go to Documents > Setup > Document Types > Diagnostic Test Result row and check Enable

5. If using Document Management > System Setup > User Groups > Document Management Tab > Patient Section > Diagnostic Test Result setting and select Full Access

6. Add Procedure Code '94762' or 'E0445' to the Item ID you will be using for the Oximetry Testing

7. Add Procedure Code "G0399" Item ID you will be using for the Home Sleep Testing; if applicable